## BD Brucella Agar with 5% Horse Blood

### USO PREVISTO

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## BIBLIOGRAFIA

1. Chu, M.C., and R.S. Weyant. 2003. *Francisella* and *Brucella*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.
2. Baron, E. J., L. R. Peterson, and S. M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.
3. Yagupsky, P. 1999. Detection of brucellae in blood cultures. J. Clin. Microbiol. 37: 3437-3442.
4. Vanderzant, C., and D. F. Splittstoesser (ed.). 1992. Compendium of methods for the microbiological examination of food, 3rd ed. American Public Health Association, Washington, D.C.
5. Hausler, W. J. (ed.). 1976. Standard methods for the examination of dairy products, 14th ed. American Public Health Association, Washington, D.C.
6. Chapin, K.C., and T.-L. Lauderdale. 2003. Reagents, stains, and media: bacteriology. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.
7. MacFaddin, J. D. 1985. Media for isolation-cultivation-identification- maintenance of medical bacteria, vol. 1, p.110-114. Williams & Wilkins, Baltimore, MD.
8. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook. American Society for Microbiology, Washington, D.C.
9. Seifert, H., et al. 1997. Sepsis – Blutkulturdiagnostik. In: MiQ - Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik, vol. 3. G. Fischer Verlag. Stuttgart, Germany.
10. Versalovic, J., and J.G. Fox. 2003. *Helicobacter*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.

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38800 Le Pont de Claix/France

Tel: +33-476 68 3636 Fax: +33-476 68 3292 <http://www.bd.com>

BD and BD logo are trademarks of Becton, Dickinson and Company.  
ATCC is a trademark of the American Type Culture Collection

Chocolate and Brucella agar w/ 5% horse blood isolates Nutritionally variant streptococci. Steak and Stab Technique of Streptococci and Enterococcus . (6 points). 1) uses 5% sheep blood agar 2) the purpose is to diagnose streptococcal pharyngitis 3) the throat swab is rolled over 1/6th of the plate to deposit the specimen 4) using loop streak inoculum. 5) stab without sterilising parts of the plate that have not been streaked, this is to create an anaerobic environment 6) the growth in the subsurface is the most reliable of haemolytic rxns (Beta) showing SLO (oxygen liable) and SLS (oxygen st .. Colony Morphology Growth seen on Blood Agar (BA) and Chocolate Agar (CA). Non-hemolytic, non-pigmented small grey colonies; may require 48-72 h for discrete colonies to become evident CO<sub>2</sub> enhances growth Slow grower; punctate colonies (0.6  $\hat{1}$ /<sub>4</sub>m - 1.5  $\hat{1}$ /<sub>4</sub>m) in 48 h. Clear or white colonies on MacConkey (MAC) agar at 48 h. Brucellosis is a commonly acquired laboratory infection; all work on suspect Brucella spp. cultures should be performed at a minimum under BSL2 conditions with BSL3 practices. Growth on BA at 48 h. Growth on CA at 48 h. . Gram Stain Tiny, often faintly stained gram-negative Blood agar is an enriched, bacterial growth medium. Fastidious organisms, such as streptococci, do not grow well on ordinary growth media. Blood agar is a type of growth medium (trypticase soy agar enriched with 5% sheep blood) that encourages the growth of bacteria, such as streptococci, that otherwise wouldn't grow. Beta hemolysis in sheep blood agar. Blood contains inhibitors for certain bacteria such as Neisseria and Haemophilus genera and the blood agar must be heated to inactivate these inhibitors and to release essential growth factors (e.g., V factor). Heating of blood agar converts it